

# AccuPen™ Dispensing Device

Our pharmacy is proud to introduce you to the patient-friendly AccuPen, an accurate dispensing device for your topical medication. The AccuPen is a unique design allowing for easier application of your measured medication.


## Directions

**Unique Locking Feature!** Turn the Push-Button End Cap clockwise to 'lock' the device. Make sure the colored arrows are not lined up.

**Simple to Use:** First, take off the removable Cap. To use the AccuPen, simply turn the Push-Button End Cap counter-clockwise to unlock the pen until the colored arrows meet, and then push the Push-Button Cap like you would push an ink pen. The medication will then be pumped out of the tip! It's that easy!

Ask your pharmacist for the AccuPen for easy and convenient dispensing of your topical medications.

Have questions? Ask the pharmacy staff – we're here to help!



Customized Medications For Your Individual Needs

138 Arnett Blvd | Sherwood Shopping Center | Danville, VA 24540

Tel: 434.792.8281 | Fax: 434.792.3235 | [pharmacist@kare-rx.com](mailto:pharmacist@kare-rx.com) | [www.kare-rx.com](http://www.kare-rx.com)